# ThermoFlash®

EC REP

**VISIONED GROUP SA**  
112, Avenue Kléber  
75116 PARIS - France

**Visiomed Technology Co., LTD.**  
2<sup>nd</sup> Floor of N°.1 building, Jia An  
Technological Industrial Park  
67 District, Bao An  
518101 SHENZHEN - China

DISTRIBUÉ PAR / DISTRIBUTED BY

Marché Européen /  
European market

**VISIONED GROUP SA**  
112, Avenue Kléber  
75116 PARIS - France  
contact@visiomed-lab.fr  
Tel : +33 (8) 92 350 334  
(0,34€/min)  
[www.visiomed-lab.com](http://www.visiomed-lab.com)

LX-260TE

By **Visiomed®**



A

2-5cm

B

Infrared Thermometer is intended for intermittent measurement of human body temperature in people of all ages without contact.

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Chamberlain, J.M., et al., Determination of Normal Ear Temperature with an Infrared Emission Detection Thermometer, *Annals of Emergency Medicine*, January 1995; / Chamberlain, J.M., et al., Determination of Normal Ear Temperature with an Infrared Emission Detection Thermometer, *Annals of Emergency Medicine*, enero de 1995; / Chamberlain, J.M., et al., Determination of Normal Ear Temperature with an Infrared Emission Detection Thermometer, *Annales de médecine d'urgence*, janvier 1995.

THERMOFLASH® LX-260TE: PRÉCISION / PRECISION / PRECISIE / PRECISIÓN / GENAUIGKEIT / PRECISIENE

34°C - 35.9°C = ± 0.3°C	93.2°F - 96.6°F = ± 0.3°F	Suivant / According to / Overeenkomstig de / Según / Gemäß / Conforme alle ASTM Standard E1965-1998 (2009)
36°C - 39°C = ± 0.2°C	96.8°F - 102.2°F = ± 0.2°F	
39°C - 42.5°C = ± 0.3°C	102.2°F - 108.5°F = ± 0.3°F	

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
# ThermoFlash®

## USER MANUAL

**THERMOFLASH® LX-260T EVOLUTION** is clinically tested and approved ASTM 1965-1998 (2009).

## CONTENTS

I. INTRODUCTION .....	17
II. PRECAUTIONS BEFORE USE .....	17
III. OPERATING PRINCIPLES.....	18
The different methods of temperature measurement .....	18
Advantages of temporal artery (TA) temperature.....	19
Practical considerations when taking a temperature .....	19
How to measure a temperature .....	19
Constraints.....	20
IV. FUNCTIONS.....	20
V. DESCRIPTION OF THE THERMOFLASH® LX-260T EVOLUTION.....	21
VI. STARTING UP YOUR THERMOFLASH® LX-260T EVOLUTION .....	21
VII. SETTINGS.....	21
VIII. TECHNICAL CHARACTERISTICS & PRECISION .....	22
IX. ADVICE .....	23
X. ACCESSORIES SUPPLIED.....	23
XI. TROUBLESHOOTING .....	23
GUARANTEE.....	69

- Follow the maintenance advice stipulated in this instruction manual.
- This device may be used for personal home use.
- This device must only be used for the purposes described in this instruction manual.
- This device must only be used in an ambient temperature range of between 10 and 40°C.
- This device must always be kept in a clean, dry area.
- Do not expose this thermometer to electric shocks.
- Do not expose this thermometer to extreme temperature conditions of >50°C or < -20°C.
- Do not use this device in relative humidity higher than 85%.
- The protective glass over the lens is the most fragile part of the thermometer.
- Do not touch the glass of the infrared lens with your fingers.
- Clean the glass with a cotton bud lightly moistened with 70° alcohol.

- Do not expose the thermometer to sunlight or to water.
- Do not use this device outside.
- Never drop the device.
- Should a problem occur with your device, please contact your retailer. Do not attempt to repair this device yourself.
- Please keep the thermometer off the children's touch. Prevent the product from children playing, swallowing and eating.
- The **MEDICAL ELECTRICAL EQUIPMENT** needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the EMC DATA.
- Portable and mobile RF communications equipment can affect **MEDICAL ELECTRICAL EQUIPMENT**.
- The measuring result of this product is only for your reference. If you have any doubt, please contact your doctor.

### Important safeguards

Check the batteries are inserted the right way.  
Reversing the polarity may damage the thermometer and affect its warranty.

*The manufacturer reserves the right to alter the specifications of the product without prior notification.*

## I. INTRODUCTION

The **ThermoFlash® LX-260T EVOLUTION**, developed by **Visiomed®**, is a revolutionary infra-red medical thermometer without contact, fitted with cutting edge MicroSecond Flash automatic calibration technology™ using latest infrared technology. This latter allows temporal artery (TA) temperature to be taken at a distance of about 2-5cm away from the forehead. Precise, instantaneous and without contact, the **ThermoFlash® LX-260T EVOLUTION** is totally suitable for a precise temperature measurement without risk. It has been shown that this TA temperature measurement, for newborns, is more precise than tympanic thermometry and better tolerated than rectal thermometry (1).

However, as with other types of thermometer, it is essential to use the **ThermoFlash® LX-260T EVOLUTION** properly in order to obtain reliable and stable results.

You are therefore advised to read this instruction manual and the safety precautions carefully before use.

(1) *Greenes D, Fleisher G. Accuracy of a Non-invasive Temporal Artery Thermometer for Use in Infants. Arch Pediatr Adolesc Med 2001; 155:376.*

## II. PRECAUTIONS BEFORE USE

The **THERMOFLASH® LX-260T EVOLUTION** is pre-set at the factory. It is not necessary to calibrate the device when starting it up.

In order to obtain reliable and stable results, you are advised each time there is a

significant change in the ambient temperature due to a change in environment, to allow the **ThermoFlash® LX-260T EVOLUTION** to acclimatise to this ambient temperature for 15 to 20 minutes before using it. It is important to allow a one minute interval between two measurements.

### III. OPERATING PRINCIPLES

All objects, solid, liquid or gas, emit energy by radiation. The intensity of this energy depends on the temperature of the object. The **ThermoFlash® LX-260T EVOLUTION** infrared thermometer is therefore able to measure the temperature of a person by the energy the person emits. This measurement can be taken thanks to an external temperature probe on the device which permanently analyses and registers the ambient temperature. Therefore, as soon as the operator holds the thermometer near the body and activates the radiation sensor, the measurement is taken instantly by detection of the infrared heat generated by the arterial blood flow. Body heat can therefore be measured without any interference from the heat of the surrounding environment.

#### The different methods of temperature measurement

##### - Core temperature

Core temperature is the most precise measurement and involves measuring the temperature in the pulmonary artery by means of a catheter equipped with a thermal probe which can read the temperature in situ. The same method is employed for probes measuring the oesophageal temperature. However, such invasive temperature measurement methods require specific equipment and expertise.

##### - Rectal thermometry

Rectal temperature adjusts slowly in comparison to the evolution of the body's internal temperature. It has been demonstrated that rectal temperature remains raised long after the internal temperature of the patient has started to drop and vice versa. Furthermore, rectal perforations have been known to occur as a result of this method and without appropriate sterilisation techniques, rectal thermometry can spread germs often found in faeces.

##### - Oral thermometry

Oral temperature is easily influenced by recent ingestion of food or drinks and by breathing through the mouth. To measure oral temperature, the mouth must remain closed and the tongue lowered for three to four minutes which is a difficult task for young children to accomplish.

##### - Axillary (armpit) temperature

Although it may be easy to measure axillary temperature, it has been proven that it does not provide an accurate measurement of the child's internal temperature. To take this type of temperature, the thermometer must be wedged tightly over the axillary artery. Despite the low sensitivity and relative inaccuracy of axillary temperature in detecting fever, this method is recommended by *The American Academy of Pediatrics* as a screening test for fever in newborns.

##### - Tympanic thermometry

In order to obtain a precise temperature reading, good command of the measurement technique is required. The thermometer probe must be placed as close as possible to the warmest part of the external ear canal. An incorrectly placed probe could lead to a false temperature reading.

##### Advantages of temporal artery (TA) temperature

Infrared arterial temperature can be measured using a device placed on the forehead, in the temporal artery region. It has been demonstrated that this relatively new method of measuring temperature is more precise than tympanic thermometry and better tolerated than rectal thermometry. The **ThermoFlash® LX-260T EVOLUTION** thermometer has been designed to produce an instant forehead temperature reading without any contact with the temporal artery. As this artery is quite close to the surface of this skin and therefore accessible and given the blood flow is permanent and regular, it allows precise measurement of the temperature. This artery is linked to the heart by the carotid artery which is directly linked to the aorta. It forms part of the main trunk of the arterial system. The efficiency, speed and comfort of taking a temperature from this area make it ideal compared with other temperature measurements methods.

##### Practical considerations when taking a temperature

- In order to ensure that precise and accurate temperature measurements are obtained, it is essential that each user has received adequate information on and training in the temperature measurement technique when using such a device.
- It is essential to remember that although procedures such as taking a temperature may be simple they must not be trivialised.
- Temperature should be taken in a neutral context. The patient must not have undertaken vigorous physical activity prior to taking his/her temperature and the room temperature must be moderate.
- Be aware of physiological variations in temperature which must be taken into consideration when evaluating the results: temperature increases by 0.5°C between 6am and 3pm. Women have a temperature that is higher, on average, by around 0.2°C. Their temperature also varies in accordance with their ovarian cycle. It rises by 0.5°C in the second half of the cycle and at the early stages of pregnancy.
- When sitting, temperature is lower by about 0.3° to 0.4°C than when standing.

##### How to measure a temperature

SEE SCHEMA  PAGE 2.

**Aim at the FOREHEAD, over the right temporal region, from a distance of about 2-5cm**, press the thermometer's measurement button and the temperature is instantly displayed.

The reliability of the measurement cannot be guaranteed if the temperature is measured over another part of the body (e.g. arm, torso...).

## Constraints

Please observe the following before any temperature measurement to ensure a stable and reliable result:

- Push back hair from the forehead.
- Wipe away any perspiration from the forehead.
- Avoid any drafts (e.g. from nasal specs, air conditioning...).
- Allow a **1 minute** interval between two measurements.
- The thermometer cannot be used as long as the display screen is not switched off.
- Each time there is a significant change in the ambient temperature due to a change in environment, to allow the **ThermoFlash® LX-260T EVOLUTION** to acclimatise to this ambient temperature for at least 15 minutes before using it.

## IV. FUNCTIONS

1. Specially designed to take the body temperature of a person regardless of the room temperature.
2. Quick and reliable results as it use the HEIMMANN infrared detection system.
3. Vocal indication of measurements and color indicators to quickly display and interpret the diagnosis of the temperature obtained.
4. Display measurement results:
  - Measurement range from 32.0 to 42.9°C (89.6 to 109.2°F).
  - Hi* displays when the result is above 42.9°C (109.2°F).
  - Lo* displays when the results is below 32.0°C (89.6°F).
5. Sound alarm if temperature is exceeded.
6. Memorization of the 32 last measures.
7. Measurement and display of room temperature:
  - Measurement ranges from 5.0°C to 59.9°C (41.0°F to 139.8°F).
8. LCD back-lighted digital screen.
9. Data displayed in Celsius or Fahrenheit.
10. Automatic stop (energy saver).
11. Small, convenient, easy to use.

### Additional uses:

The **ThermoFlash® LX-260T EVOLUTION** can also be used to measure the temperature of a baby-bottle or bath, or room temperature (by using the *SURFACE TEMP* or *ROOM* function). This function is in accordance with the Directive 89/336/EEC Electromagnetic Compatibility.

## V. DESCRIPTION OF THE THERMOFLASH® LX-260T EVOLUTION

SEE SCHEMA **B** PAGE 2.

- |                         |                           |                          |
|-------------------------|---------------------------|--------------------------|
| <b>1</b> Infra-red lens | <b>2</b> LCD screen       | <b>3</b> Measure button  |
| <b>4</b> Battery case   | <b>5</b> Infra-red sensor | <b>6</b> Setting buttons |

## VI. STARTING UP YOUR THERMOFLASH® LX-260T EVOLUTION

1. Install battery
2. For the first use or when inserting new batteries wait between 10 minutes for the warm-up of the apparatus and when inserting the new batteries.
3. Aim towards the forehead (see **A** for the **ThermoFlash® LX-260T EVOLUTION** positioning), from a distance of 2-5 cm, press the measuring key, the temperature is displayed in 1 second. The temperature can also be taken behind the ear lobe.
4. Before taking the temperature, make sure to remove hair and perspiration from the forehead. You can also take the temperature behind the ear lobe.
5. Replace new battery when you see the lower battery logo (  ).

## VII. SETTINGS

### 1. Choosing the temperature mode

The **ThermoFlash® LX-260T EVOLUTION** is specially designed to measure body temperature without contact, of children or adults. It can be used by consumers in household environment and doctors in clinic as reference. For this, use the *BODY* mode.

Measurement range for *BODY* mode: 32°C – 42.9°C (90°F – 109°F).

*Note : You can also use the ThermoFlash® LX-260T EVOLUTION to measure the temperature of an area or an object, food, a liquid or a room temperature. For this, use the SURFACE TEMP mode.*

Measurement range for *SURFACE TEMP* mode: 0°C – 60°C (32°F – 140°F).

By selecting the *ROOM* Mode, it will display the ambient temperature of the room you are in. Press the *MODE* button for one second, the screen will display: *BODY* for the body measurement temperature. Press the *MODE* button again to switch to *SURFACE TEMP* or *ROOM* mode.

*Note: The ThermoFlash® LX-260T EVOLUTION is automatically set to BODY.*

**Important:** The surface temperature differs from the internal body temperature. To obtain the internal temperature, always use the *BODY* mode.

Please make sure to select the *BODY* mode for an internal temperature reading and the *SURFACE TEMP* mode or the *ROOM* mode for an external area reading (bottle, bath, room...).

### 2. Choosing the temperature unit

Press the °C/°F button to select the temperature unit. The symbol of the chosen temperature unit is displayed (°C for Celsius degrees and °F for Fahrenheit degrees).

### 3. Data memory

To display the last temperature measurement, press the *MEM* button.

You will obtain the last temperature measurement. To switch on the last measurement, press the °C/°F button.

The number indicated in interval of two measures corresponds to the number of measurements. Press on the button *MODE* for review of the previous measurements.

To exit data memory, press the measure key. The **Thermoflash® LX-260T EVOLUTION** switches off after 5 seconds.

### 4. Choosing the language

Use the Language key to choose a language for the vocal function. The choice of language is visible.

Languages available: French, English, Spanish, Italian, German, Dutch, Silent Function.

By default the **ThermoFlash® LX-260T EVOLUTION** is in French.

### 5. Alarm setup

The alarm threshold default value is 38°C (100.4°F). It cannot be changed. In case of higher temperatures than 38°C, a signal beeps while measuring the temperature.

### 6. Changing the batteries

Display: when the LCD screen displays Battery, the battery is used.

Operation: Open the lid and change the batteries, taking great care with the correct positioning. A mistake with this could cause damage to the apparatus and compromise the guarantee of your **ThermoFlash® LX-260T EVOLUTION**. Never use rechargeable batteries. Use only batteries for single usage.

## VIII. TECHNICAL CHARACTERISTICS & PRECISION

1. Normal conditions of use operating temperature: 10°C ~ 40°C (50°~ 104°F). Normal conditions of storage and transportation temperature: -20°C ~ 50°C (68°~ 122°F). Humidity rate: 10 - 85%
2. Power: DC3 V (2 batteries LR6).
3. Size: 220 x 145 x 48 mm.
4. Weight: 180g.
5. Display Resolution: 0.1° C (0.1°F).
6. Measuring range:
  - In *BODY* mode: 32°C – 42.9°C (90°F – 109°F).
  - In *SURFACE TEMP* mode: 0°C – 60°C (32°F to 140°F).
7. Precision: From 36°C to 39°C (96.8°F to 102.2°F) = +/-0.2°C.
8. Consumption: ≤ 50mW.
9. Measuring distance: 2-5cm.
10. Automatic stop: 5 sec.

### Longevity use

The **ThermoFlash® LX-260T EVOLUTION** is guaranteed for 40.000 readings.

### IX. ADVICE

- The protective glass over the lens is the most important and fragile part of the thermometer, please take great care of it.
- Clean the glass with cotton fabric, wet with water or 70° alcohol.
- Do not use other batteries than mentioned batteries, do not recharge non rechargeable batteries, do not throw in fire.
- Remove the batteries when thermometer is not used for an extended period of time.
- Do not expose the thermometer to sunlight or water.
- An impact will damage the product.

### X. ACCESSORIES SUPPLIED

User manual, Guarantee card, Batteries supplied (LR6).

### XI. TROUBLESHOOTING

If you have one of the following problems while using your **ThermoFlash® LX-260T EVOLUTION**, please refer to this breakdown service guide to help resolve the problem. If the problem persists, please contact our customer service at +33 892 350 334.

- **The screen displays temperature superior to 95°:** the temperature is in Fahrenheit. Change the measurement unit to Celsius.
- **The screen displays the body temperature inferior to 32°C (89.6°F):** to measure a body temperature, the *BODY* mode must be selected. If you are on *SURFACE TEMP* mode, the 32°C (89.6°F) temperature displayed is showing the external temperature that your body releases.
- **The screen displays the message HI:** when using the **ThermoFlash® LX-260T EVOLUTION**, the message *HI* can show on the screen. The analysis is above the measurement range selected, either superior to 42.9°C (109°F) in *BODY* Mode or superior to 60°C (140°F) in *SURFACE TEMP* Mode.
- **The screen displays the message LO:** when using the **ThermoFlash® LX-260T EVOLUTION**, the message *Lo* can show on the screen. The temperature analysed is under the measuring range selected, either less than 32°C (90°F) in *BODY* Mode or less than 0°C (32°F) in *SURFACE TEMP* Mode.

This message displays in various cases – please find below a list of the main cases:

Reasons for LO message display	Advice
Temperature reading hampered by hair, perspiration...	Make sure that there is no obstruction prior to taking a temperature.
Temperature hampered by an air flux.	Make sure there is no air flux as this could interfere with the infrared system.
Temperature readings too close together, the ThermoFlash® LX-260T EVOLUTION did not have the chance to boot itself.	Respect the pause of 15 seconds minimum between two readings – 1 minute pause is advised.
The measuring distance is too far.	Please respect the measuring distance (between 2 and 5 cm).

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EMC DATA

## EMC data

**Guidance and manufacturer's declaration - electromagnetic immunity - for equipment and system that are not life-supporting.** The model LX-260TE is intended for use in the electromagnetic environment specified below. The customer or the user of the model LX-260TE should assure that it is used in such an environment.

Immunity test	EN 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF	3 Vrms	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the model LX-260TE, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[ \frac{3.5}{E1} \right] \sqrt{P}$ 80 MHz to 800 MHz  $d = \left[ \frac{7}{E1} \right] \sqrt{P}$ 800 MHz to 2.5 GHz  where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in meters (m). Fields strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
IEC 61000-4-6	150 kHz to 80 MHz		
Radiated RF	3 V/m	3 V/m	
IEC 61000-4-3	80 MHz to 2.5 GHz		

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.

a. Fields strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model LX-260TE is used exceeds the applicable RF compliance level above. The model LX-260TE should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary such as reorienting or relocating the model LX-260TE.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

**Recommended separation distances between portable and mobile RF communications equipment and the equipment or system - for equipment and systems that are not life-supporting.** The model LX- 260TE is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the model LX-260TE can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the model LX-260TE as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output of transmitter	Separation distance according to frequency of transmitter m	
	80 MHz to 800 MHz $d = \left[ \frac{3.5}{E1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[ \frac{7}{E1} \right] \sqrt{P}$
0.01	0.12	0.23
0.1	0.38	0.73
1	1.2	2.3
10	3.8	7.3
100	12	23

For transmitters rated at a maximum output power not listed above the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagations is affected by absorption and reflection from structures, objects and people.

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**EN:** Visiomed\* will repair or replace this product free of charge in the case of defective parts or manufacturing defects, in accordance with the conditions mentioned below as follows:

DURATION: 24 MONTHS RETURN TO WORKSHOP

**LIMITS AND EXCLUSIONS:** This guarantee concerns only the original final purchaser. A purchase invoice, or another proof of purchase, with this guarantee card will be required to obtain an after-sales service, in accordance with this guarantee. This guarantee card will not be extended to another person only the original final purchaser. This guarantee becomes void if the serial numbers on the product are modified, replaced, illegible, absent, or if repair has been carried out by a service not approved, including the user.

This guarantee covers only the defects of the material or parts, occurring during normal use of the product. It does not cover the damage caused during the transport of the apparatus, causes due to repairs being carried out by the distributor, by any modifications undertaken, any connection of equipment not approved by Visiomed\*, or causes contrary to those written in the user manual or notice. Moreover, the present guarantee does not cover damage due to falls, bad handling, bad installations, damage by fire, floods, lightning, or any other natural disaster. This guarantee does not cover the packing of the material, the accessories, the defects caused by commercial exposure of the product, show room, sale space, demonstration etc... Normal maintenance, cleaning and the replacement of parts where wear is normal, are not covered by the terms of this guarantee. Visiomed\* and its representatives and agents will not in any case be held responsible for any damage and consecutive damages due to the mishandling of this product. This guarantee is the only valid one at Visiomed\*, any other guarantee (commercial guarantee) except this one will not be taken into account.

**IMPORTANT:** During the guarantee period if you are dissatisfied with the repairs of this product, please contact the Visiomed\* customer service.

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**EN:** *Disposal of Waste Equipment by Users in Private Households in the European Union.*

*The symbol on the product or on its packaging indicates that this product must not be disposed of with your other household waste. Instead, it is your responsibility to dispose of your waste equipment by handing it over to a designated collection point for the recycling of waste electrical and electronic equipment. The separate collection and recycling of your waste equipment at the time of disposal will help to conserve natural resources and ensure that it is recycled in a manner that protects human health and the environment. For more information about where you can drop off your waste equipment for recycling, please contact your local city office or your household waste disposal service.*

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© 03/2017 LX-260TE
**EN:** *TensioFlash is a brand registered by Visiomed Group SA. The constructor keeps the right to modify the photographs and characteristics without previous notice. Registered design. The logos and trademarks are the property of their respective owners.*